



**Department of Pathology  
Point of Care  
Urine Drug Screen Reporting Form**

Test Performed by:  _____ (Name)	Date of Service, Ordering Provider Patient Name <small>Attach Patient Demographics Label</small> Patient MRN, DOB <small>(at least 2 patient identifiers)</small> Patient Sex, Age, CSN
_____	_____
(Name)	(Time)

Adulterants						
<b>Creatinine</b> (circle one) (Read at 45 seconds)	<b>Abnormal</b>		<b>Normal</b>			
	Neg	10	20	50	100	200
<b>Specific Gravity</b> (circle one) (Read at 45 seconds)	<b>Abnormal</b>		<b>Normal</b>			
	1	≥1.030	1.005	1.015	1.025	
<b>Oxidant/Bleach</b> (circle one) (Read at 30 seconds)	Positive			Negative		
<b>Drug</b> (Read at 5 minutes)	<b>Result</b> (POS/NEG)			<b>Control</b> (✓)		
<b>AMP</b> Amphetamines *Cut-off Level 1000 ng/ml						
<b>BAR</b> Barbiturates *Cut-off Level 300 ng/ml						
<b>BUP</b> Buprenorphine *Cut-off Level 10 ng/ml						
<b>BZO</b> Benzodiazepine *Cut-off Level 300 ng/ml						
<b>COC</b> Cocaine *Cut-off Level 300 ng/ml						
<b>MET</b> Methamphetamine *Cut-off Level 1000 ng/ml						
<b>OPI300</b> Morphine *Cut-off Level 300 ng/ml						
<b>OXY</b> Oxycodone *Cut-off Level 100 ng/ml						
<b>MTD</b> Methadone *Cut-off Level 300 ng/ml						
<b>THC</b> Marijuana *Cut-off Level 50 ng/ml						

Fentanyl is not detected in this assay.

Time in hours since last dose: \_\_\_\_\_

Record results in the correct box above.

**Positive** - enter POS. Result is positive when there is **ONE** line is seen (control line) and no band in

**Negative**- enter NEG. Result is negative when there are **TWO** lines seen (control and test line).

**Control line present, enter ✓ .**

If no control line present the test is invalid, enter **INVALID**. Repeat test with new iCup.