

*University of Michigan Health System  
Molecular Diagnostics Laboratory*

**Microsatellite Instability (MSI) Testing by PCR Slide Process Form**

Name \_\_\_\_\_ Reg # \_\_\_\_\_

Block # \_\_\_\_\_ Specimen site \_\_\_\_\_

Date received \_\_\_\_\_

---

---

**Histology Sectioning Request for MSI:**

Date received: \_\_\_\_\_

**Sectioning Instructions:**

10 micron sections

5 unstained slides + 1 H&E per block

**DO NOT BAKE THE SLIDES!**

**Once complete, please take all slides, block, and paperwork to 2F331D (GA signout room).**

Processed by \_\_\_\_\_ Date: \_\_\_\_\_


---

---

**GA Fellow/ Resident Slide Review**

Please clearly mark on the H&E slide the normal tissue area (solid line) and tumor tissue area (dotted line) as shown below. There should be sufficient delineation between each. The normal tissue can be tested from a separate block from the same patient if necessary.

**Once complete, take all slides and paperwork to Specimen Processing to be forwarded to Molecular Diagnostics.**



Completed by \_\_\_\_\_ Date: \_\_\_\_\_